Infection Control Division | **Accessories**

## IC10/20FR Auto-Reader.

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators and Hygiene Monitoring Systems.

### Usage

For incubation and readout of:

Bionova® Rapid Fluorescence Readout Biological Indicators: BT220 (Steam): 3 h at 60 °C, BT221 (Steam): 1 h at 60 °C and BT110 (EO): 4 h at 37 °C, BT102 (Formaldehyde): 2 h at 60 °C, BT95 (VH2O2): 2 h at 60 °C.

Bionova® Super Rapid Fluorescence Readout Biological Indicators: BT222 (Steam): 1 h at 60 °C, BT223 (Steam): 30 min. at 60 °C, BT96 (VH2O2): 30 min. at 60 °C.

Bionova® Ultra Rapid Fluorescence Readout Biological Indicators: BT224 (Steam): 20 min. at 60 °C.

Chemdye® Protein Pen: 15 min. at 60 °C

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Low Voltage Directive 2014/35/EU.

Electromagnetic Compatibility Directive 2014/30/EU.

RoHS Directive 2011/65/EU.

WEEE Directive 2012/19/EU.

### FDA 510(k)

K163646

### Characteristics

12 metal positions (3.8 cm deep, 0.9 cm diameter) to incubate Rapid/Super Rapid/Ultra Rapid Fluorescence Readout Self-contained biological indicators.

1 metal position (3.8 cm deep, 1.0 cm diameter) to incubate Protein Pen.

Special hole for external thermometer.

1 ampoule crusher.

Dual Temperature System: Allows selecting between two different incubation temperatures.

Allows running different incubation times simultaneously.

USB connection for registering results in a PC through Reading and Traceability Software.

Thermal Printer to record results in sterilization book.

Automatic detection and cancelation of biological indicator fluorescence reading.

Fluorescence reading wavelength: ~ 460 nm

Dimensions: 18.0 cm high, 26.0 cm diameter.

Voltage range: 100 – 240 V AC.

Power: 28 W.

Frequency: 50-60 Hz.

### Environmental conditions during manufacture

T = 15-30 °C, RH = 35-85 %.

Infection Control Division | **Accessories**

## IC10/20FR Auto-Reader.

For incubation and readout of Rapid, Super Rapid or Ultra Rapid Biological Indicators and Hygiene Monitoring Systems.

**Storage conditions**

T = 10-30 °C, RH = 30-80 %, keep away from light.

**Transport conditions**

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

**Shelf life**

5 years.

**Packing**

1 unit per box.

Packing information: product code and description, storage conditions and manufacturer information.

**Labelling**

On the product and on product's box: product code and description, storage conditions, batch number, serial number, barcode and manufacturer information.

**Possible target markets**

Healthcare, Food, Pharmaceutical and medical devices industries.

**Other important information**

IC10/20FR has a button located on the incubator panel to select temperatures at 37 °C or 60 °C. It also has a button which allows setting an incubation time for half of the positions (1 to 6), and another button to select a different incubation time for the remaining positions (7 to 12).

IC10/20FR Incubator has an automatic readout system and a printer to record the results.

**Precautions**

Do not pour any liquid inside.

Do not immerse into any liquid.

Only for indoor use.

Disconnect the power cord before cleaning.

Do not use abrasive or corrosive cleaners or disinfectants.

Make sure that the incubator is connected to a properly rated power cord.

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