

Sterilisation Wraps



Our single-use sterilisation wraps help to prevent microorganisms (bacteria, viruses, spores) from contaminating tray contents. They provide effective sterilisation, provide physical protection and maintain sterility until the point of use and support.

Our technologies:

- Crepe paper
- SMS
- Wetlaid non-wovens



Application:

- Sterilisation wraps for terminally sterilised medical devices
- Container filters for effective sterilisation protection
- Overwrap for single-use Custom Procedure Trays (CPT)
- Tray liner designed to absorb moisture and protect trays
- Sterilisation wraps for terminally sterilised medical devices
- Filter for sterilisation containers
- Overwrap for single-use Custom Procedure Trays (CPT)
- Tray liner providing moisture absorption and tray protection

Sterilisation Methods:

The Sterile Barrier System products are designed for wrapping medical devices that can be terminally sterilised using the following methods:



STEAM



EO



FORM



IRRAD



H2O2

✉ hello@sychem.co.uk

☎ 02380 240729

🌐 sychem.co.uk