



SYCHEM CHANGE

Our class-leading dual access change station

Sychem provides a comprehensive range of animal care solutions, services & consumables to support research facilities in providing the best conditions for the animals in their care.



Active Airflow Technology

SychemCHANGE utilises Sychem's Active Airflow Technology to protect both animals and technicians.

A filtered laminar airflow creates a downflow of H14 HEPA filtered air onto the worksurface.

ISO Class 3 airflow ensures a barrier of protection between animal and operator, limiting the potential for cross-contamination and maximising staff safety.

Designed as a dual-access workstation, the unit is compatible with all rodent cages and features a fully reversible TRESPA® worksurface, available in either a fixed raised or lowered position.



Easy Cage Change

SychemCHANGE has a large work area, ideal for various rodent cage sizes. The cage lid holder supplied as standard enables quick, aseptic cage changing. The station is perfect for complete, base-only cage changing and cleaning procedures.



Ergonomic Excellence

SychemCHANGE offers the best height range on the market, making it suitable for any user. At the touch of a button, the electrically driven legs can be raised or lowered, maximising the user's comfort levels and simplifying the cleaning process.



TRESPA® worksurface

SyhemCHANGE features a fully reversible TRESPA® worksurface, providing a chemical & impact resistant worktop. Clients can order the worksurface in either a fixed raised or lowered position.



SYCHEMCHANGE



AIRFLOW RANGE

SychemCHANGE is part of our Airflow Range, a trio of systems that also includes SychemFILL – a free flow bedding dispenser and SychemWASTE – a bedding disposal unit.







