




CLEANING INDICATORS

At Sychem, in partnership with **Terragene**, we offer a complete portfolio of Cleaning Indicators.

- CDWA4
- Cleaning Indicators
- Partnership
- Popular Cleaning Indicators
- Indicator holders

 Contact us today!
0845 644 6824

 hello@sychem.co.uk  sychem.co.uk

Sychem and Terragene have worked in partnership since 2017.



CDWA4

Cleaning Indicator

Our high-level performance Indicator for cleaning procedures. The CDWA4 can detect all factors that affect the result of a washing process.

- Challenge level: very high
- Allows detection of all factors affecting the outcome of cleaning processes
- Initial Colour: fuchsia
- Size: 70 mm x 30 mm
- Synthetic temperature-stable support
- Result: elimination of fuchsia organic substance
- To interpret results visually, use the results reference guide on product's instructions for use
- Optionally, the readout of the cleaning indicator can be done using Terragene's Trazanto Automatic Quality Control and Traceability System

Cleaning Indicators


Cleaning Indicators are used for validation and routine monitoring of efficacy of cleaning processes in washers, particularly in the Health industry.

- Washing Indicators
- Endoscope Process Challenge devices
- Thermodisinfection Indicators





UK Exclusive partnership with Terragene



As the UK exclusive partner for Terragene since 2017, Sychem provides an array of infection prevention and control solutions for hospitals, dental clinics, food preparation facilities, the cosmetic industry, pharmaceutical sites and medical devices.

Some Cleaning Indicators can be used in conjunction with the Trazanto- a unique AI device for quality and traceability control of washing, disinfection, and sterilisation processes.



Popular Cleaning Indicators

CDWA3 Splat: A Cleaning Indicator for cleaning performance tests.

Lumenia SixFlow: Challenge kit monitoring cleaning of endoscopes in Automatic Endoscope Reprocessors.

IT27W-5: Thermodisinfection Indicator for moist heat disinfection.

CDWU: A cleaning indicator for monitoring ultrasonic washing machines.



Indicator holders

CDWAH: A holder for Washing Indicators

CDWAH-U: A holder for Ultrasonic Washing Indicators



Chemdye®
Lumenia



Contact us today!

0845 644 6824